



Hands on Experiment :  
Determining the **enthalpy change**  
of a reaction with the use of  
Arduino system

# Setup for the Experiment


Determining the enthalpy change of metal displacement reaction



Using Arduino system with temperature sensor to record the temperature change during the reaction



# Equipment required for Experiment




Zinc powder



Copper (II) sulphate



# Equipment required for Arduino system



Waterproof DS18B20  
Temperature Sensor


3 Jumper wires (Male to Male)




Resistor (4.7k ohm)

Arduino UNO

Solderless Breadboard

# Basic : Solderless breadboard(麵包板)






| Positive      | Negative       |
|---------------|----------------|
| Power         | Ground         |
| Plus sign (+) | Minus sign (-) |
| Red           | Blue or black  |




# Operation of the Arduino system

## 1. Wiring of temperature sensor in Arduino system




# Identify the pins

Signal : 2




Power : 5V (+)  
GND (-)

Connect red jumper wire to 5V




Connect black jumper wire to GND




Connect Yellow jumper wire to 2



Connect another end of red jumper wire to breadboard +




Connect another end of black jumper wire to breadboard -




Bend the 4.7k ohm resistor to “C” shape




Connect the banded resistor to breadboard + and A12




Connect another end of yellow jumper wire to breadboard B12




Connect red and black wire of temp sensor to breadboard + and - respectively



Connect yellow wire of temp sensor to breadboard E12



Check all wires and resistor to see all are firmly secured on Arduino and breadboard



<https://www.youtube.com/watch?v=p6AN6rPWzfl>

# Experimental Setup

Copper(II)sulphate solution  
Zinc powder

Clamp and stand


Arduino System with  
temperature sensor

Polystyrene cup


Magnetic Stir bar

Magnetic Stirrer

Computer



# Points to note



Avoid glass rod from contact with temp sensor



Magnetic Stir bar should not to contact with temp sensor

## 2. Copy the following into Arduino Libraries



One Wire

Dallas Temperature\_Library

Waterproof DS18B20  
Temperature Sensor



Program Files (x86)>Arduino>libraries



# 3. Upload Arduino code



Waterproof DS18B20  
Temperature Sensor

Open Arduino IDE software

Copy the code into Arduino IDE

Save the Arduino Sketch as “temp”



```
temp | Arduino 1.8.1
File Edit Sketch Tools Help
temp
#include <OneWire.h>
#include <DallasTemperature.h>

// Arduino數位腳位2接到1-Wire裝置
#define ONE_WIRE_BUS 2

// 應用程式庫建立物件
OneWire oneWire(ONE_WIRE_BUS);
DallasTemperature sensors(&oneWire);

void setup(void) {
  Serial.begin(9600);
  Serial.println("Temperature Sensor");
  // 初始化
  sensors.begin();
  Serial.println("CLEARDATA");
  Serial.println("LABEL,Time,Temperature");
}

void loop(void) {
  Serial.print("DATA,TIME,");



  // 要求匯流排上的所有感測器進行溫度轉換
  sensors.requestTemperatures();

  // 取得溫度讀數 (攝氏) 並輸出，
  // 參數0代表匯流排上第0個1-Wire裝置
  Serial.println(sensors.getTempCByIndex(0));


  //delay(1000);
}
```

Arduino/Genuino Uno on COM8

# 4. Checking the port of Arduino system to record the temperature data



# 5. Uploading of code into Arduino Board



```
temperature | Arduino 1.8.1
File Edit Sketch Tools Help
temperature
#include <OneWire.h>
#include <DallasTemperature.h>

// Arduino數位腳位2接到1-Wire裝置
#define ONE_WIRE_BUS 2

// 運用程式庫建立物件
OneWire oneWire(ONE_WIRE_BUS);
DallasTemperature sensors(&oneWire);

void setup(void) {
  Serial.begin(9600);
  Serial.println("Temperature Sensor");
  // 初始化
  sensors.begin();
  Serial.println("CLEARDATA");
  Serial.println("LABEL,Time,Temperature");
}

void loop(void) {
  Serial.print("DATA,TIME,");

  // 要求匯流排上的所有感測器進行溫度轉換
  sensors.requestTemperatures();

  // 取得溫度讀數 (攝氏) 並輸出,
  // 參數0代表匯流排上第0個1-Wire裝置
}
```

Arduino/Genuino Uno on COM4



```
temperature | Arduino 1.8.1
File Edit Sketch Tools Help
temperature
#include <OneWire.h>
#include <DallasTemperature.h>

// Arduino數位腳位2接到1-Wire裝置
#define ONE_WIRE_BUS 2

// 運用程式庫建立物件
OneWire oneWire(ONE_WIRE_BUS);
DallasTemperature sensors(&oneWire);

void setup(void) {
  Serial.begin(9600);
  Serial.println("Temperature Sensor");
  // 初始化
  sensors.begin();
  Serial.println("CLEARDATA");
  Serial.println("LABEL,Time,Temperature");
}

void loop(void) {
  Serial.print("DATA,TIME,");

  // 要求匯流排上的所有感測器進行溫度轉換
  sensors.requestTemperatures();



  // 取得溫度讀數 (攝氏) 並輸出,
  // 參數0代表匯流排上第0個1-Wire裝置
}
```

Done uploading.

Sketch uses 5566 bytes (17%) of program storage space. Maximum is 32256 bytes.  
Global variables use 299 bytes (14%) of dynamic memory, leaving 1749 bytes for local variables. Maximum

Arduino/Genuino Uno on COM4


# 6. Checking the port of Arduino system to start record the temperature data by PLX-DAQ



|    | A       | B           | C |
|----|---------|-------------|---|
| 1  | Time    | Temperature |   |
| 2  | 9:11:51 | 24.06       |   |
| 3  | 9:11:51 | 24.06       |   |
| 4  | 9:11:51 | 24.06       |   |
| 5  | 9:11:51 | 24.06       |   |
| 6  | 9:11:51 | 24.06       |   |
| 7  | 9:11:51 | 24.06       |   |
| 8  | 9:11:52 | 24.06       |   |
| 9  | 9:11:52 | 24.06       |   |
| 10 | 9:11:52 | 24.06       |   |
| 11 | 9:11:52 | 24.06       |   |
| 12 | 9:11:52 | 24.00       |   |
| 13 | 9:11:52 | 24.00       |   |
| 14 | 9:11:52 | 24.00       |   |
| 15 | 9:11:53 | 24.06       |   |
| 16 | 9:11:54 | 24.00       |   |
| 17 | 9:11:54 | 24.06       |   |

# 7.Data Collection by Excel plug-ins

|    | A           | B                  |
|----|-------------|--------------------|
| 1  | <b>Time</b> | <b>Temperature</b> |
| 2  | 上午 10:11:43 | 20.19              |
| 3  | 上午 10:11:45 | 20.19              |
| 4  | 上午 10:11:47 | 20.19              |
| 5  | 上午 10:11:48 | 20.19              |
| 6  | 上午 10:11:50 | 20.19              |
| 7  | 上午 10:11:52 | 20.19              |
| 8  | 上午 10:11:54 | 20.19              |
| 9  | 上午 10:11:56 | 20.19              |
| 10 | 上午 10:11:57 | 20.19              |
| 11 | 上午 10:11:59 | 20.19              |
| 12 | 上午 10:12:01 | 20.19              |
| 13 | 上午 10:12:03 | 20.19              |
| 14 | 上午 10:12:05 | 20.19              |
| 15 | 上午 10:12:06 | 20.19              |
| 16 | 上午 10:12:08 | 20.19              |
| 17 | 上午 10:12:10 | 20.19              |
| 18 | 上午 10:12:12 | 20.19              |
| 19 | 上午 10:12:13 | 20.19              |
| 20 | 上午 10:12:15 | 20.19              |
| 21 | 上午 10:12:17 | 20.19              |
| 22 | 上午 10:12:19 | 20.19              |
| 23 | 上午 10:12:21 | 20.19              |
| 24 | 上午 10:12:22 | 20.19              |
| 25 | 上午 10:12:24 | 20.19              |



Total duration : 540 seconds    540 data

# 8. Data analysis by using Microsoft Excel

|    | A           | B         | C              | D     | E         | F          | G                     | H                    | I  | J | K | L   | M |
|----|-------------|-----------|----------------|-------|-----------|------------|-----------------------|----------------------|--|---|---|---|---|
| 1  | Actual Time | Time diff | Duration (sec) | Temp  | Initial T | Decrease T |                       |                      |  |   |   | Temperature Data Collected by Arduino UNO |   |
| 2  | 10:11:43    | 0:00:00   | 0              | 20.19 | 20.19     |            |                       |                      |  |   |   | Date of Expt: 13/2/2017                   |   |
| 3  | 10:11:45    | 0:00:02   | 2              | 20.19 | 20.19     |            |                       |                      |  |   |   | Sample Rate : 1 second                    |   |
| 4  | 10:11:47    | 0:00:04   | 4              | 20.19 | 20.19     |            |                       |                      |  |   |   | Reaction Enthalpy                         |   |
| 5  | 10:11:48    | 0:00:05   | 5              | 20.19 | 20.19     |            |                       |                      |  |   |   | Zn(s) + CuSO4(aq) → ZnSO4(aq) + Cu(s)     |   |
| 6  | 10:11:50    | 0:00:07   | 7              | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 7  | 10:11:52    | 0:00:09   | 9              | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 8  | 10:11:54    | 0:00:11   | 11             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 9  | 10:11:56    | 0:00:13   | 13             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 10 | 10:11:57    | 0:00:14   | 14             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 11 | 10:11:59    | 0:00:16   | 16             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 12 | 10:12:01    | 0:00:18   | 18             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 13 | 10:12:03    | 0:00:20   | 20             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 14 | 10:12:05    | 0:00:22   | 22             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 15 | 10:12:06    | 0:00:23   | 23             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 16 | 10:12:08    | 0:00:25   | 25             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 17 | 10:12:10    | 0:00:27   | 27             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 18 | 10:12:12    | 0:00:29   | 29             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 19 | 10:12:13    | 0:00:30   | 30             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 20 | 10:12:15    | 0:00:32   | 32             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 21 | 10:12:17    | 0:00:34   | 34             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 22 | 10:12:19    | 0:00:36   | 36             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 23 | 10:12:21    | 0:00:38   | 38             | 20.19 | 20.19     |            |                       |                      |  |   |   |   |   |
| 24 | 10:12:22    | 0:00:39   | 39             | 20.19 | 20.19     |            | Data Collection :     | Calculations :       |  |   |   |   |   |
| 25 | 10:12:24    | 0:00:41   | 41             | 20.19 | 20.19     |            | Metal :               | Zinc                 | mass of solution (m) = 50 g  |   |   |   |   |
| 26 | 10:12:26    | 0:00:43   | 43             | 20.19 | 20.19     |            | Solution :            | Copper (II) sulphate | Sp Heat Cap of solution (c) = 4.18 J g <sup>-1</sup> K <sup>-1</sup> |   |   |   |   |
| 27 | 10:12:28    | 0:00:45   | 45             | 20.19 | 20.19     |            | Time for adding zinc: | 173 Sec              | Temp Diff (ΔT) 48.87 °C  |   |   |   |   |
| 28 | 10:12:29    | 0:00:46   | 46             | 20.19 | 20.19     |            | Maximum Temp          |                      | heat capacity of container (c') 0 J K <sup>-1</sup>                  |   |   |   |   |
| 29 | 10:12:31    | 0:00:48   | 48             | 20.19 | 20.19     |            | y =                   | -0.0469 x + 77.171   |  |   |   |   |   |
| 30 | 10:12:33    | 0:00:50   | 50             | 20.19 | 20.19     |            | T2 =                  | 69.06                |  |   |   |   |   |
| 31 | 10:12:35    | 0:00:52   | 52             | 20.19 | 20.19     |            |                       |                      | Number of mole   |   |   |   |   |
| 32 | 10:12:37    | 0:00:54   | 54             | 20.19 | 20.19     |            | Initial Temp          |                      | zinc 0.0616 mol  |   |   |   |   |
| 33 | 10:12:38    | 0:00:55   | 55             | 20.19 | 20.19     |            | T1 =                  | 20.19                | Copper (II) sulphate 0.0500 mol                                      |   |   |   |   |
| 34 | 10:12:40    | 0:00:57   | 57             | 20.19 | 20.19     |            |                       |                      | Limiting number of mole 0.0500 mol                                   |   |   |   |   |

# Scientific Inquiry opportunity

To investigate the enthalpy change of neutralization